



ITO QCM SENSOR CHIP

Diameter 14mm (for QCM-I and other instruments)



FRONT side (sensing side)



REAR side (electrical contact side)

| | |
|---|----------------------------|
| Nominal frequency: | 4950±50KHz |
| Crystal cutting: | AT |
| Material: | Pure Z quartz |
| Angle: | 35° 15'±30' |
| XX' Angle: | 0°00'±30' |
| Quartz diameter: | 14mm |
| Electrode on the sensor side (wrap around): | ITO |
| Electrode on the contact side: | ITO (see below) |
| Appearance : | Transparent in the visible |